

An easy guide to making your ZomaJet® Pen growth hormone (10mg)

Before you're ready to Jet, Set, Grow! you'll need to make your growth hormone. Every week you'll need to get a new batch ready.

Put a new adaptor on the Zomacton® (10mg) medicine bottle

- Remove the cap
- Wipe the rubber seal with a sterile wipe
- On a flat surface, PUSH the adaptor down



Make your growth hormone

- Remove the rubber cap from the pre-filled syringe, and take the cap off the adaptor
- PUSH the syringe into the top of the adaptor, and gently pull the plunger up a bit to free any air bubbles
- PUSH DOWN the plunger slowly until it hits the bottom
- PULL the syringe away. Some of the adapter will come away with it
- Put the protective cap back onto the adaptor
- Your growth hormone is now ready to use!



Keep the medicine bottle nice and cold in the fridge

It can stay fresh for 28 days

Your ZomaJet® Pen

